


Cerrahi KS-W Serisi, çok parlak optik lens benimseyerek aşırı parlaklık, üstün kaliteli düzgün ışınlar sunan taşınabilir LED Cerrahi far sistemlerinin en popüler nesillerinden biridir.

KS-W02 black



Dual Lithium Polymer battery packs provide constant steady brightness and longevity to the headlight. A full charge illuminates for up to 5+5 hours at the highest brightness

Built-in Rechargeable Lithium Polymer battery with LED screen display power indication, fast and easy recharging of battery

Precision multiple lens optical system consisting of four high-quality lens. Completely homogeneous illumination of the whole field of view, and a clearly defined light spot. Spot size adjustable by inner mini light grating



Optional yellow filter can effectively reduce blue light and avoid light curing in dental surgery.



100-240V wide range compatible transformer

LED Medical Headlight KS-W series		
Model	KS-W01	KS-W02
Light source	Single LED	
LED Power	3W	5W
LED span life	M 5000hours	
Illumination	18000-41000lu	20000-48000lu
Work Distance	^ 200-500m	^ 200-500m
Light spot	Variable from D=20mm to D=110mm	
Power Supply	AC 90-240V / DC 5V	
Battery NO.	Single(Standard),Dual(optional)	
Capacity of battery	1800mAh per piece	
Charging Time	2 hours	
Battery using time	5 hours(max illumination)-20hours(min illumination)	
Brightness control	Digital button on/off and brightness control with LED	
Color of appreance	Black or White	
Optional Magnifier	M250 2.5X (Black&White)/ M350 3.5X (Black&White)/ M450 2.5-3.5X Black / M-4X,K-4X 4X Black/ K-5X 5X Black/ K-6X 6X Black	
Package	EPE Foam package with color box (standard)/ Alluminium case (option)	



EPE Foam inside



Color box (Standard)



Options: Loupes, Battery, Alluminium case